

# **ELISA Buffer Set (Phosphate system,96T)**

Catalog Number: EBS-002

Pack Size: 96 tests

IMPORTANT: Please carefully read this manual before performing your experiment.

<u>For Research Use Only. Not For Use in Diagnostic or Therapeutic Procedures</u>



# **MATERIALS PROVIDED**

Table 1. Materials provided

Catalog	Components	Size (96T)	Format	Storage	
				Unopened	Opened
EBS002-C01	High-bind plate	1 plate	Powder	2-8°C	2-8°C
EBS002-C02	Coating buffer	12 mL	Liquid	2-8°C	2-8°C
EBS002-C03	10xWashing Buffer	60 mL	Liquid	2-8°C	2-8°C
EBS002-C04	Blocking buffer	50 mL	Liquid	2-8°C	2-8°C
EBS002-C05	Substrate solution	12 mL	Liquid	2-8°C, avoid light	2-8°C, avoid light
EBS002-C06	Stop Solution	7 mL	Liquid	2-8°C	2-8°C

# **SRORAGE**

- 1. The unopened kit is stable for 12 months from the date of manufacture if stored at 2°C to 8°C.
- 2. The opened kit should be stored per Table 1. The shelf life is 30 days from the date of opening.

Note: a. Do not use reagents past their expiration date.

- b. Find the expiration date on the outside packaging.
- c. The reagent should be properly balanced to room temperature before use.

## **APPLICATION**

**Asia and Pacific:** 

- 1. High-bind plate (EBS002-C01): 96 well microplates. Used for solid phase adsorption of samples in coated solution.
- 2. Coating buffer (EBS002-C02): Coated buffer for diluting coated sample.
- 3. 10xWashing Buffer (EBS002-C03): For washing the enzyme label plate, 10xWashing Buffer should be diluted with purified water to 1xWashing Buffer before use (For example, 60 mL 10xWashing Buffer+ 540 mL purified water).
- 4. Blocking buffer (EBS002-C04): Used to block the enzyme label plate. Diluting the Blocking buffer 4 times with 1xWashing Buffer can be used as a Sample buffer diluent (For example, 1 mL Blocking buffer + 3 mL 1xWashing Buffer).

1 / 4

US and Canada: Tel: +1 800-810-0816

Tel: +86 400-682-2521

E-mail: order@acrobiosystems.com

Web: http://www.acrobiosystems.com



ACTO\*

5. Substrate solution (EBS002-C05): For color reaction.

6. Stop Solution (EBS002-C06): Used to terminate the experiment.

# **TYPICAL DATA**

## 1. Coating

Dilute Human TNF-α Capture Antibody stock solution (ACROBiosystems Cat.No.CRS-D002, CRD002-C01) to 1.0 μg/mL with Coating Buffer (EBS002-C02) to make Human TNF-α Capture Antibody working solution.

Description 100 μL of Human TNF-α Capture Antibody working liquid (1.0 μg/mL) is added to each well onto a High-bind plate (EBS002-C01), sealed with a microporous plate film, and incubated at 4°C overnight (or for 16 hours).

## 2. Washing

Remove the remaining solution by aspiration, add 300 µL of 1×Washing Buffer to each well, gently tap the plate for 1 minute, remove any remaining 1×Washing Buffer by aspirating or decanting, invert the plate and blot it against paper towels. Repeat the wash step above for three times.

#### 3. Blocking

Add 300  $\mu$ L Blocking Buffer (EBS002-C04) to each well, seal the plate with microplate sealing film and incubate at room temperature for 2.0 hours.

## 4. Washing

Repeat step 2.

## 5. Add Samples

The reconstructed Human TNF- $\alpha$  Standard (ACROBiosystems Cat.No.CRS-D002, CRD002-C02) was diluted with Sample buffer, and the dilution range was 4.89-312.5 pg/mL. Add 100  $\mu$ L Samples to each well. For blank Control wells, please add 100  $\mu$ L Sample buffer.

#### 6. Incubation

Seal the plate with microplate sealing film and incubate at room temperature for 1 hour.

#### 7. Washing

Repeat step 2.

2 / 4

US and Canada:

Tel: +1 800-810-0816

Web: http://www.acrobiosystems.com



## 8. Add Human TNF-α Detection Antibody

Dilute Biotinylated-Human TNF- $\alpha$  Detection Antibody (ACROBiosystems Cat.No.CRS-D002, CRD002-C03) stock solution to 0.5  $\mu$ g/mL with Sample buffer to make Biotinylated-Human TNF- $\alpha$  Detection Antibody working solution.

For all wells, add 100  $\mu$ L Biotinylated-Human TNF- $\alpha$  Detection Antibody (0.5  $\mu$ g/mL) working solution. Please prepare it for one-time use only.

## 9. Incubation

Seal the plate with microplate sealing film and incubate at room temperature for 1 hour.

## 10. Washing

Repeat step 2.

# 11. Add Streptavidin-HRP

For all wells, add 100 μL Streptavidin-HRP (ACROBiosystems Cat.No.CRS-D002, CRD002-C04) (dilute at 1:2000) working solution. Please prepare it for one-time use only, avoid light.

#### 12. Incubation

Seal the plate with microplate sealing film and incubate at room temperature for 30 min.

## 13. Washing

Repeat step 2.

#### 14. Substrate Reaction

Add 100 µL Substrate Solution (EBS002-C05) to each well. Seal the plate with microplate sealing film and incubate at room temperature for 20 min, avoid light.

#### 15. Termination

Add 50 µL Stop Solution (EBS002-C06) to each well, and tap the plate gently to allow thorough mixing. *Note: The color in the wells should change from blue to yellow.* 

#### 16. Data Recording

Read the absorbance at 450 nm and 630 nm using UV/Vis microplate spectrophotometer within 10 minutes.

Note: To reduce the background noise, subtract the value read at OD450nm with the value read at OD630 nm.


3 / 4

 US and Canada:
 Tel: +1 800-810-0816
 Web: <a href="http://www.acrobiosystems.com">http://www.acrobiosystems.com</a>

 Asia and Pacific:
 Tel: +86 400-682-2521
 E-mail: <a href="mailto:order@acrobiosystems.com">order@acrobiosystems.com</a>



# 17. Data Analysis



Immobilized Human TNF-α Capture Antibody (Cat.No.CRS-D002, CRD002-C01) at 1 μg/mL (100 μL/well) can bind Human TNF-α Standard, and then add Biotinylated-Human TNF-α Detection Antibody (Cat.No.CRS-D002, CRD002-C03) at 0.5 µg/mL (100 µL/well). Detection was performed using HRP-conjugated streptavidin with a linear range of 4.89-312.5 pg/mL.

Web: http://www.acrobiosystems.com