

# Anti-Clumping Additive

**Cat. # AC-1112-11**

## For Research Use Only

Anti-Clumping Additive is a culture supplement that facilitates suspension by significantly reducing aggregation and adhesion of cells. Cell densities can achieve a higher level in optimized culture with Anti-Clumping Additive. This product is a gentle, animal origin-free, protein-free, serum-free and chemically defined formulation that minimize the interference and risk in cell cultivation. Please notice that Anti-Clumping Additive can only be used after cationic lipid transfection process for this category of cell line.

### Product Features:

- Reducing aggregation and adhesion of cells
- Facilitating suspension
- Achieving a higher cell densities
- Animal Origin-Free, Protein-Free, Serum-Free and Chemically Defined

### Storage:

Store Anti-Clumping Additive at 2°C - 8°C in the dark.

### General Specification:

Category	Parameters
Form	Liquid
Volumn	100ml
Cell Line:	HEK293, CHO, Insect cells
Shipping Condition:	Room Temperature

For detailed information or for further support, please contact our technical support team at:

**[TechSupport@acrobiosystems.com](mailto:TechSupport@acrobiosystems.com)**