

Biotinylated Monoclonal Anti-Mouse IL-10 Antibody, Rat IgG2b (AS544)

Catalog # FABr017-B



Application

Flow Cytometry.

Specificity

This product is a specific antibody specifically reacts with mouse IL-10 protein.

Immunogen

Purified mouse IL-10 Protein.

Conjugate

Biotin

Formulation

Supplied as 0.2 μ m filtered solution in PBS, pH7.4 with BSA as protectant.

Contact us for customized product form or formulation.

Storage

This product is stable after storage at:

- Store at 2° to 8°C.
- Do not use after the expiration date shown on the label.

Background

Detection antibody for Multiplex Bead Assay.

Clinical and Translational Updates

Discounts, Gifts,
and more!

